

IMProv Video Tutorial

by

<https://www.msstudio.ca/>

dschriem@ucalgary.ca
<https://www.msstudio.ca/mss-improv/>

FADE IN - msstudio IMProv project

GETTING STARTED

NARRATOR

In order to prepare the IMProv deployment bundle and execute a modelling run. We are going to look at a sample project available from github to get things started.

PREPARE IMP TOPOLOGY AND CONFIG FILES

NARRATOR

Pull up the MassSpecStudio IMProv project online tutorial

https://github.com/pellst/imp_msstudio_init/blob/master/IMPov_msstudio_tut.md

NARRATOR

This presents us with the steps that need to be taken in order to pull together the various data files (fasta, pdb, em, xl etc.). The wizard enables us to define the Topology file representation of, for example, the protein structures involved. We also need a configuration file to define the various parameters, such as the number of cpu's for the MPI job run, the number of frames as well as names and characteristics of the structures being included.

https://raw.githubusercontent.com/pellst/imp_msstudio_init/master/uml_activity_diag_improv.svg

The Activity Diagram for Integrative Modeling using IMP presents a diagrammatic overview of the Integrative Modeling lifecycle and the Stages involved.

NARRATOR

Stage 1 is what we are covering here. This involves defining the various data files that need to be included. This includes preparing the YAML config file and Topology file.

Stage 2 and 3 cover the inner workings of a IMP modeling run and the generation of the rmf files for viewing with Chimera.

Stage 4 is the Analysis phase.

NEW INTEGRATIVE MODELLING PLATFORM PROJECT

NARRATOR

Start MassSpecStudio. Open a new project and select the Integrative Modelling Platform Project template. Give your project a name (eg: IMP_PRC2) and optionally choose an alternative location.

ADD PROTEINS

NARRATOR

Using the Add Proteins wizard screen. Select the FASTA or PDB files to add reference sequences. This will then show the Name and give the opportunity to customize the Topology by clicking the Manage button under the Topology column for the row with a protein name.

ADD PROTEIN TOPOLOGY

NARRATOR

The Topology record can be edited to set the start and end of the sequence together with the PDB Offset etc. Once you click Ok you will be returned to the Add Proteins wizard screen so that you can do the same for each of the Proteins involved. The representation can be adjusted e.g. two structures can be assigned to a single sequence and bead size can be adjusted. Once you have completed all the Proteins that you wish to amend. You can click the Next button (at the bottom right hand corner of the screen) which will take you to the Add Link Data wizard screen

ADD LINK DATA

NARRATOR

Add Link Data wizard screen is where you can add additional data files including Cross-Linking, Hydrogen Exchange, Covalent Labeling and Electron Microscopy. These files will be included in their respective folders for the final output that is generated. Once you have completed your file selections you can click the Next button (at the bottom right hand

(MORE)

NARRATOR (cont'd)
corner of the screen). This will
take you to the Configure IMP
wizard screen.

HX-XL CLASSIFICATION

NARRATOR
Using the histogram view we are
able to adjust the range suitable
for setting the distance restraints
aided by the HX information. This
is where we set the bin sizes to
capture the five categories
covering: Very Loose, Loose,
Medium, Tight, Very Tight.

CONFIGURE IMP

NARRATOR
The Configure IMP wizard screen is
where we define the Directory path
to export the data files and
modeling scripts to. We also set
the Sampling Frames and States
here. The Ridgid Body and Super
Ridgid Body assignments are
available through the pick lists
provided. The final step is to
click the Export button (at the
bottom right hand corner of the
screen). This will produce the
folder structure containing the
Topology and YAML Config file
together with the raw data files
that you selected in the wizard
steps (data folder). It also adds
a folder with the modeling scripts
needed (imp_model) to perform the
job run using the python driver
script provided.

AMENDMENTS

NARRATOR
An existing project can be reopened
so that we can make adjustments to
the parameters and settings. This
enables us to tweek the Topology
and YAML configuration files using
the wizard steps again. For
advanced users the Topology and
YAML configuration files can be
opened in a text editor and
manually changed. These manual
adjustments would not be reflected
back in the original project and
may be overwritten in the event
(MORE)

NARRATOR (cont'd)
 that amendments are later made
 through the IMProv project wizard.
 Take care to keep a separate folder
 with your manual changes to prevent
 loss thereof.

Note

FADE OUT - msstudio IMProv project

Note

FADE IN - IMProv deployment

DEPLOYMENT STEPS

NARRATOR

Pull up the IMProv script
 deployment, online tutorial. We
 will start with deploying to Cedar
 on the Compute Canada cluster. AWS
 deployment will be covered later.
 In both cases the pre-requisite
 steps to setup an account and login
 to those services is mentioned in
 the online tutorial and includes
 links to their getting started
 guides.

https://github.com/pellst/imp_msstudio_init/blob/master/IMProv_uuml_diag.png

DEPLOYMENT ON CEDAR

https://github.com/pellst/imp_msstudio_init/blob/master/IMProv_on_Cedar_tut.md

NARRATOR

We make use of a setup script from
 github gist in order to provide the
 commands needed to get the sample
 project from github. This brings
 with it the example files and
 scripts that we will be using to
 complete this demonstrating.

https://gist.githubusercontent.com/pellst/4853822ea5ca74785af61d0ad39cf84d/raw/uoc_mss_prep_step1.sh

https://github.com/pellst/imp_msstudio_init/blob/master/mss_out/imp_model/uoc_mss_prep_step1.sh

DEPLOYMENT ON AWS

https://github.com/pellst/imp_msstudio_init/blob/master/IMProv_on_AWS_tut.md

NARRATOR

Note

FADE OUT - IMProv deployment

GLOSSARY

ABBREVIATIONS

Cryo-EM| cryoelectron microscopy |
<https://www.sciencedirect.com/science/article/pii/S0304416517302374>

FDR| False Discovery Rate |
<https://www.bioinfor.com/fdr-tutorial/>

HPC| High Performance Computing |
https://docs.computecanada.ca/wiki/Getting_started

HX-MS| Hydrogen eXchange Mass Spectrometry |
[https://neu.hxms.com/research/tutorial_theory.htm#:~:text=Hydrogen%20exchange%20\(HX\)%20combined%20with,of%20proteins%20and%20protein%20structure.](https://neu.hxms.com/research/tutorial_theory.htm#:~:text=Hydrogen%20exchange%20(HX)%20combined%20with,of%20proteins%20and%20protein%20structure.)

IMP| Integrative Modeling Platform |
<https://integrativemodeling.org/>

PMI| Python Modeling Interface |
<https://integrativemodeling.org/>

PRC2| Polycomb Repressive Complex 2 |
<https://www.ncbi.nlm.nih.gov/pmc/articles/PMC5008062/>

SLURM| Simple Linux Utility for Resource Management |
https://en.wikipedia.org/wiki/Slurm_Workload_Manager

XL-MS| Crosslinking Mass Spectrometry |
<https://www.technologynetworks.com/proteomics/articles/cross-linking-mass-spectrometry-a-key-player-in-the-structural-biologists-toolbox-322446>

FASTA| The FASTA format is sometimes also referred to as the "Pearson" format (after the author of the FASTA program and ditto format). |
https://www.bioinformatics.nl/tools/crab_fasta.html |
https://en.wikipedia.org/wiki/FASTA_format

PDB| The Protein Data Bank (pdb) file format is a textual file format describing the three-dimensional structures of molecules held in the Protein Data Bank |
<https://pdbe101.rcsb.org/learn/guide-to-understanding-pdb-data/introduction> | <https://www.rcsb.org/>

AWS| Amazon Web Services |
<https://aws.amazon.com/console/> | Cedar| Compute Canada HPC Cluster | <https://status.computecanada.ca/>

Linux| Operating System, RedHat Enterprise Linux (or variants, such as CentOS or Scientific Linux)